



Points to Practice

DISINFECT

- Disinfect the silicone septum with an appropriate antiseptic agent.
- Disinfect prior to each access.



ACCESS



Syringe:

- Access the BD Q-Syte device with an ISO luer lock or luer slip syringe.
- Connect using a straight-on approach.



Administration Set:

- When using an administration set with a spin nut, insert the slip luer into the BD Q-Syte device. Then secure the connection with the spin nut.



Note: Do not use needles or blunt cannula with the BD Q-Syte device.



Blood Withdrawal:

- Use ISO syringe or BD Luer-Lok™ access device.

FLUSH



- Flush according to facility policy or until device is clear. If unable to clear replacement is recommended.
- To minimize reflux of blood into the vascular access device, it is recommended to clamp before disconnecting.
- When disconnecting, hold the BD Q-Syte device firmly in place to prevent accidental removal.



The BD Q-Syte device has a unique luer penetrating septum design that provides higher flow rates, fluid path visualization, and consistent dependable performance.

NOTE: Refer to package insert for precautions and complete instructions for use.



Helping all people
live healthy lives

BD Medical
Sandy, Utah
84070
888.237.2762
www.bd.com