


PRODUCT DATA SHEET

PROSET INTRAFIX SAFESET(4186097)



Attribute Name	Attribute Value
Product line	IV
Type of set	IV admin. set
Intended use	for infusion therapy
Bestseller	No
IV-Set: Type	Intrafix SafeSet
IV-Set: Spike	integrated
IV-Set: Application	Gravity
IV-Set: Filter	Fluid filter 15 µm
IV-Set: AirStop	Yes
IV-Set: PrimeStop	Yes
IV-Set: Add. Accesses	-
IV-Set: Number Add. Accesses	0

Attribute Name	Attribute Value
IV- Set tube: length	180 cm
IV-Set tube: Inner-ø	3.00 mm
IV-Set tube: Outer-ø	4.10 mm
IV-Set tube: Flow control	Roller clamp
Access 1 Inline: BCV	No
Access 1 Inline	No
Access 1 Inline: Number	0
Access 1 Inline: Type	-
Access 1 Inline: Accessories	-
Extension: Inner-ø tube mm	0.00 mm
Extension: Outer-ø tube mm	0.00 mm
Extension: length	0 cm
Extension: Type	-
Type of filter	No filter
Access 2 patient: BCV	Yes
Access 2 patient	Yes
Access 2 patient: Number	2
Access 2 patient: Type	Discifix
Access 2 patient: Accessories	-
Distal Tube: Length	10 cm
Distal Tube: Inner-ø	2.00 mm
Distal Tube: Outer-ø tube	4.10 mm
Fixed connection	Alcohol connection
Length total	190 cm
length category	150 - 300 cm
Material (Tube)	PVC
Light-protected	No
Stock-managed	Yes

Attribute Name	Attribute Value
Market limitation	No
Customer limitation	No
Sterility	sterile
Disposable / Reusable	disposable
Product contains Latex	No

