





Product data sheet - 06.23.12

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.48 (Plastic), 450.00.54 (OMF Instr), 450.00.57 (Plastic), 450.00.58 (Cranialis), 450.00.59 (Spinalis), 450.00.60 (Plastic neu), <b>Brochures</b> 451.02.03 (Gubisch), 451.04.76 (ARIS), <b>Single sheets</b> 451.06.85 (Brustdiss),		
Last update	2008-01-24		
Title	<b>Delicate tissue forceps</b>		
Article number	<b>06.23.12</b>	Author	Micro-Adson
Description	Delicate tissue forceps		
Length	12 cm - 4 3/4	Weight	
		Diameter	
Properties	straight	Tip / Working end	1 x 2 teeth
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826007041		
GMDN	11785		
Risk category EC	1		
UMDNS	14-257		
Net weight	0.0200		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 06.23.12

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4024	X15Cr13	C 0,12 - 0,17% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---