





III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 07.82.07

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1),		
Last update	2008-01-24		
Title	Splinter Forceps		
Article number	07.82.07	Author	Feilchenfeld
Description	Splinter Forceps		
Length	8 cm - 3 1/8	Weight	
		Diameter	
Properties	straight, serrated	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826010041		
GMDN	11777		
Risk category EC	1		
UMDNS	14-257		
Net weight	0.0100		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 07.82.07

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4024	X15Cr13	C 0,12 - 0,17% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---