





Product data sheet - 88.12.06

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.58 (Cranialis), 450.00.59 (Spinalis),		
Last update	2013-09-23		
Title	Coagulation Forceps, bipolar, straight, blunt		
Article number	88.12.06	Author	
Description	Coagulation Forceps, bipolar, straight, blunt		
Length	11 cm - 4 3/8"	Weight	
		Diameter	
Properties	sterelizable insulating plastic coating	Tip / Working end	0.6 mm
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	G880702_2014_12		
GTIN	4046826163259		
GMDN	47848		
Risk category EC	2b		
UMDNS	11-502		
Net weight	0.0200		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 88.12.06

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material	RILSAN		

Delivery mode	
---------------	---