

Product Specifications

Brand : Klinion® Easy Care
Product : Vaginal Speculum

Packaging & Product Reference

Product Reference	Description	Piece(s) per Box	Box Dimension (LxWxH in cm)*
531400	KEC VAGINAL SPEC 20MM XS WHITE	200	36 x 80 x 39,5
531401	KEC VAGINAL SPEC 24MM S BLUE	200	36 x 80 x 39,5
531402	KEC VAGINAL SPEC 26MM M RED	200	36 x 80 x 39,5
531403	KEC VAGINAL SPEC 30MM L GREEN	200	36 x 80 x 39,5

* Sizes are approximate

Product description

Klinion Easy Care disposable Sterile vaginal speculum.

Intended use

The Klinion Easy Care vaginal speculums are intended to be used for single use vaginal examination.

Classification

Class IS under Council Directive MDD 93/42/EEC

Features


- Available four sizes distinguishable through colour: white, blue, red, green
- Smooth surface and rounded edges
- Single use
- Manufactured with body temperature friendly material
- Allows the adjustment of the vaginal opening

Product Specifications

Product composition

- Polystyrene - Colour : transparent
- This product is:
 - Latex-free
 - Phthalates-free

Specifications

Drawing	REF Color	Dimension (mm)			
		531400 WHITE	531401 BLUE	531402 RED	531403 GREEN
	Tip Width (W)	20	24	26	30
	Penetration Length	106	106	106	106
	Maximum Length (L)	115	115	115	115

Sterilization

Sterilized with ethylene oxide.

Storage conditions


Store in a dry place. Do not expose to the direct sunlight or to sources of heat.

Shelf life

5 years.

Product Specifications

Picture




Ref. 531400

Ref. 531403

Ref. 531402

Ref. 531401

Approval	Name	Date	Signature
Product Manager	Levine van Buren	13 September 2018	
RA Manager	Francis Dekker	13 September 2018	