




# Texium<sup>®</sup> Syringe

Reducing the hazards of dangerous drugs



## No risk of disconnections

- Eliminates any risk of a syringe accidentally disconnecting from the Texium<sup>®</sup> Closed Male Luer.
- The Texium Syringe consists of a **syringe permanently connected** to the Texium closed male luer device and it is designed to enable the safe handling of fluids and medications, particularly hazardous or cytotoxic drugs.

## No spinning = Ease of handling

- Holding the Texium syringe to inject or withdraw creates a more robust system that makes it easier and quicker to use (for example with over-pressurised bags or collapsible container).
- There is no need to hold the Texium closed male luer when connecting/disconnecting from the needle free valve or female luer.
- Texium Syringe makes it easier to withdraw the desired amount of fluid.

## No need for adapters

- No cumbersome and **costly** extra adaptors are needed
- It can be connected to any ISO compliant female luer.
- Texium will ensure a closed system when connected to a SmartSite<sup>®</sup> needle free valve.

## Sealing the system for hazardous drugs

- Needles and vials present a high risk of leakage. The Texium Syringe with the SmartSite<sup>®</sup> needle-free valve delivers leak-free connections and disconnections. This significantly reduces surface contamination and accidental exposure to hazardous drugs.

Reconstitution >>> Preparation >>> Administration

Chemo safety system



# What you can't see can hurt you



## Intended use

The Texium Syringe is designed to provide leak-free handling of potentially hazardous fluids, such as chemotherapy, radioactive isotopes, and blood products, during all three phases of drug handling when the Texium Syringe is used in conjunction with the SmartSite Needle Free Valve port: the preparation of the drug, the administration of the drug to the patient, and waste handling.

## Ordering Information

Model Number	Description	Case quantity
MY8003-0006	Texium 3 ml Needle-Free Syringe	100
MY8005-0006	Texium 5 ml Needle-Free Syringe	100
MY8010-0006	Texium 10 ml Needle-Free Syringe	100
MY8020-0006	Texium 20 ml Needle-Free Syringe	100
MY8060-0006	Texium 60 ml Needle-Free Syringe	50

Easier stock management with one model code

[carefusion.com](http://carefusion.com)

